

MEDICAL INSIGHTS

Bronchoscopic removal of blood clots in the central airway

Role of cryoextraction for restoring of airway patency

Background

The COVID-19 pandemic has increased the global number of patients requiring intensive care including mechanical ventilation – sometimes with ECMO support. Blood clot formation with life-threatening airway obstruction has been reported in critically-ill patients on ECMO¹. Other causes include tuberculosis, bronchiectasis, mycetoma, cancer, cystic fibrosis and bleeding after interventions (e.g. biopsies)^{3,8,9}.

Although there is no data on the incidence of airway obstruction by blood clots in the patient population severely affected by COVID-19 as yet, blood clot removal is a topic of high interest at the interface between intensive care medicine and bronchoscopy.

Method

A reasearch of the currently available literature was conducted by Erbe using PubMed, Livivo and Google Scholar. The terms "blood clot AND bronchial", "lung AND cryo" and "blood clot AND cryo AND broncho" were searched.

Challenges and goals

Blood clot removal can be achieved with different methodologies. However, guidelines as to when to choose which treatment option are not yet present.

Flexible bronchoscopy is now the main method used and can be carried out at bedside^{3,8}. Most blood clots can be removed with suction, lavage and flexible forceps³. However, these methods include the risk of bleeding and airway trauma and might not be successful due to the fragile structure of blood clots¹¹⁻¹³.

Rigid bronchoscopy requires general anesthesia and bedside application can be challenging^{2,3,14,15}. The use of balloon catheters increases the risk of bronchial injury and mucosa damage, and topical thrombolytic agents pose a risk for rebleeding^{5,14}.

With cryoextraction, large casts can be removed and bedside application is possible². One case report describes bleeding after cryoextraction¹⁶. Current literature confirms a rapid removal of blood clots with flexible cryoprobes where other methods have failed³.

Nevertheless, it is currently perceived as a back-up solution.

Results and key findings

Cryoextraction represents a standard treatment for blood clot removal in the airway according to the guidelines^{6,7}. The success rate of cryoextraction ranges between 87.5% and 100%¹. En bloc removal was frequently achieved. Extraction in a piecemeal fashion was also described¹ and the overall procedure time was between 7 and 15 minutes^{3,5}.

Conclusion


Cryoextraction of blood clots is a safe procedure^{1,3}. Advantages of cryoextraction include performance with either flexible or rigid intubation. Absence of general anesthesia in flexible bronchoscopy can be beneficial in critically-ill patients, which makes cryoextraction an option in several procedural environments, including bedside application. Friable blood clots were described to be stabilized by freezing².

Taking into account current data with regards to safety, speed, recanalization rate and procedural environment, cryoextraction can be considered as a *primus inter pares* (first among equals) of the available options for blood clot removal.


Products

The evaluated studies were conducted with ERBECRYO[®] 2 and ERBOKRYO CA units together with flexible cryoprobes. To support the standardization of blood clot removal, single use cryoprobes with 1.7 mm (20402-410) and 2.4 mm (20402-411) can be used¹⁰.

Freezing times between 10 and 20 seconds were reported^{3,4}.



CRYO.ERBE-MED.COM



References

- Schmidt LH et al. Blood clot removal by cryoextraction in critically ill patients with pulmonary hemorrhage. *J Thorac Dis*. 2019 Oct;11(10):4319-4327. doi: 10.21037/jtd.2019.09.46.
- Sriratanaviriyakul N et al. Safety and Clinical Utility of Flexible Bronchoscopic Cryoextraction in Patients With Non-neoplasm Tracheobronchial Obstruction: A Retrospective Chart Review. *J Bronchology Interv Pulmonol*. 2015 Oct;22(4):288-93. doi: 10.1097/LBR.000000000000203.
- Lee H et al. Successful removal of endobronchial blood clots using bronchoscopic cryotherapy at bedside in the intensive care unit. *Tuberc Respir Dis (Seoul)*. 2014 Oct;77(4):193-6. doi: 10.4046/trd.2014.77.4.193. Epub 2014 Oct 31.
- Cook MI et al. The use of cryotherapy via bronchoscopy for removal of obstructing tracheobronchial thrombi. *Int J Crit Illn Inj Sci*. 2015 Jul-Sep;5(3):215-6. doi: 10.4103/2229-5151.165009.
- Sehgal IS et al. Use of a Flexible Cryoprobe for Removal of Tracheobronchial Blood Clots. *Respir Care*. 2015 Jul;60(7):e128-31. doi: 10.4187/respcare.03861.
- Du Rand IA et al. British Thoracic Society Interventional Bronchoscopy Guideline Group. British Thoracic Society guideline for advanced diagnostic and therapeutic flexible bronchoscopy in adults. *Thorax*. 2011 Nov;66 Suppl 3:iii1-21. doi: 10.1136/thoraxjnl-2011-200713.
- Bolliger CT et al. ERS/ATS statement on interventional pulmonology. European Respiratory Society/American Thoracic Society. *Eur Respir J*. 2002 Feb;19(2):356-73. doi: 10.1183/09031936.02.00204602.
- Davidson K, et al. Managing Massive Hemoptysis. *Chest* 2020; 157: 77–88
- Radchenko C, Alraiyes AH, Shojaee S. A systematic approach to the management of massive hemoptysis. *J Thorac Dis* 2017; 9: S1069-S1086
- Hetzel J et al. Evaluation of Efficacy of a New Cryoprobe for Transbronchial Cryobiopsy: A Randomized, Controlled in vivo Animal Study. *Respiration*. 2020;99(3):248-256.
- Dobbertin I, Dierkesmann R, Kohlhäufel M. *Lehrbuch und Atlas der Bronchoskopie: Kapitel 17. Hämoptoe, pathologische Gefäßbefunde*. 2nd ed.; 2008
- Bhardwaj B, Bhardwaj H, Youness H et al. Bronchoscopic cryoextraction: a novel approach for the removal of massive endobronchial blood clots causing acute airway obstruction. *Southwest J Pulm Crit Care* 2013; 7: 184–189
- Engelhardt K, Pirolli T, Raman L et al. Successful Use of Pulmonary Cryotherapy for Tracheobronchial Thrombus Extraction and Recanalization of the Tracheobronchial Tree During a Pediatric Venovenous Extracorporeal Membrane Oxygenation Run. *Pediatr Allergy Immunol Pulmonol* 2019; 32: 28–30
- Saladi L, Lvovsky D. Organized blood clot masquerading as endobronchial tumor: A review of management and recent advances. *Respir Med Case Rep* 2018; 24: 165–169
- Özdemir C, Sökücü SN, Kocatürk Cİ et al. Are flexible bronchoscopic cryoextraction practices an alternative to rigid bronchoscopy? *Tuberc Toraks* 2019; 67: 15–21
- Tendă ED, Yakub A, Pitoyo CW et al. Combination of bronchoscopic cryoextraction and argon plasma coagulation in treatment of total central airway obstruction caused by giant blood clot formation in massive airway bleeding. *Respir Med Case Rep* 2016; 19: 9–11

Important information

We have prepared this document with care. Nonetheless, we cannot completely rule out errors in this document.

The information, recommendations and other data ("Information") contained in this document reflect our state of knowledge and the state of science and technology at the time of preparing the document. The information is of a general nature, non-binding and serves solely for general information purposes and does not represent instructions for use or notes on application.

The information and recommendations contained in this document do not constitute any legal obligations on Erbe Elektromedizin GmbH or their associated companies ("Erbe") or any other claims against Erbe. The information does not represent a guarantee or other quality statement; these require an express contractual arrangement with Erbe in individual cases.

Erbe shall not be liable for any type of damage resulting from following information given in this document, regardless of the legal reason for liability.

Every user of an Erbe product is responsible for checking the respective Erbe product for its properties as well as the suitability for the intended type of application or intended purpose in advance. The suitable type of application of the respective Erbe product is given by the user manual and the notes on use for the corresponding Erbe product. The user is obliged to check whether the existing user manual and the notes on use correspond with the status for the specific Erbe product. The devices may only be used according to the user manual and the notes on use.

The information on setting values, application sites, duration of application and the use of the respective Erbe product is based on the clinical experience of physicians independent from Erbe. They represent guidelines which need to be checked by the user for their suitability for the actual planned application. Depending on the circumstances of an actual application case, it may be necessary to deviate from the information provided. The user is responsible for checking this in each case when using an Erbe product. We wish to point out that science and technology is constantly subject to new developments arising from research and clinical experience. For this reason, it may be necessary for the user to deviate from the information provided in this document.

This document contains information about Erbe products which may possibly not be approved in a specific country. The user of the respective Erbe product is obliged to inform him/herself as to whether the Erbe product he/she is using is legally approved in his/her country and/or if legal requirements or restrictions for use possibly exist and to what extent.

This document is not intended for users in the USA.

Erbe Elektromedizin GmbH
Waldhoernlestrasse 17
72072 Tuebingen
Germany

Phone +49 7071 755-0
info@erbe-med.com
erbe-med.com
medical-videos.com